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1. General

Please read this document carefully before using the system and keep it for future consultation! Failure to comply with these instructions for use may cause damage to patient's health, implants and instruments.

This document does not contain all of the information necessary for the selection and application of the system. For safe and proper handling, refer to further product-related instructions, such as the surgical technique associated with the system, the packaging and cleaning instructions and the device labels on the packaging. For instruments that require calibration, please refer to the Additional Product Information, Torque Wrench (IFU-EU-134-140-001).

2. User Group and Environment

The products may only be used and operated in an aseptic medical environment by persons who have the appropriate training, knowledge or experience in the orthopedic and surgical field.

The intended users are physician, theatre nurse as well as CSSD staff (CSSD = Central Sterile Supply Department).

3. Patient Group

The patient group for our medical device is made up of adult, anesthetized patients of any ethnic origin and of any gender, in whom the described indications are present and who are not considered unsuitable on the grounds of the listed contraindications. The device is made for older, less mobile or immobile patients.

4. System Description

The product is part of a mechanical reconstruction of the hip joint.

The Modular Trauma Heads are large metal heads that are part of a hemi hip arthroplasty often performed in emergency care cases after femoral neck fractures. They are used in combination with a femoral stem. The system consists of two parts, the heads that articulate in the native acetabulum and sleeves that built the conjunction between the heads and the hip stem.

The instruments are designed to define the implant size of the Modular Trauma Heads. The set consists out of the trial heads and the neck length adapter. There are four different neck length adapters for the Modular Trauma Heads. To determine the size of the head there is a handle available, which can be connected to the trial heads.

The extractor for the sleeves can be used when a sleeve has to be removed from a head or stem.

The existing instruments enable the surgeon a size measurement of the head size as well as a trial reduction with an adjustment of the neck length.

5. Intended Use

The non-active, surgically-invasive implantable Modular Trauma Head manufactured by Waldemar Link GmbH & Co. KG is intended for long-term replacement of the femoral side of a diseased and / or defective hip joint in the human body. It replaces the prosthesis head of the human hip. The Modular Trauma Head forms a hemiarthroplasty of the hip joint when combined with the prosthesis stem. The Modular Trauma Head can be used with full-grown, anesthetized patients of any ethnic origin and sex. The Modular Trauma Head is implanted without cement.

The implants may only be used and operated in an aseptic medical environment by persons who have the required training, knowledge and experience in the orthopedic and surgical field. The implants are supplied in sterile condition individually packed as single-use products.

The purpose of the instrumentation is to allow the user to use the associated implant system as part of the procedures described in the associated surgical technique. Any other use of the instrumentation is not permitted.

The instruments consist of defined, combinable instruments.

All instruments of the instrument set are intended for transient use.

6. Indications

General Indications Implants:

- Mobility-limiting diseases, fractures or defects of the hip joint or proximal femur which cannot be treated by conservative or osteosynthetic procedures.

Indications:

- Femoral neck fractures ¹⁾

¹⁾ For older, less mobile or immobile patients

General Indications Instruments:

- The Instruments for the Modular Trauma Heads in the present combination are exclusively intended for the application and implantation of the Modular Trauma Head.

7. Contraindications

Implants:

- Acute and chronic infections, local and systemic insofar as they compromise the successful implantation of a hemiarthroplasty
- Allergies to (implant) materials
- Insufficient / inadequate bone mass or quality, which prevents a stable anchorage of the prosthesis (eg. Osteoporosis, protrusio acetabuli)
- Acetabular injuries (e.g. fractures)
- Acetabular arthritis

Instruments:

- The Instruments for the Modular Trauma Heads in the present combination are not intended for the application and implantation of another medical device aside from the Modular Trauma Heads.
- Material intolerance to the instrument materials
- Instruments are not allowed to be implanted

8. Possible Risks and Side Effects

- Infection
- Instability, Dislocation
- Periprosthetic fracture
- Residual complaints
- Wear

9. Clinical Benefit

The clinical benefit of the Modular Trauma Heads is defined as:

- Restoration of joint functionality.
- Reduced time of surgery, compared to conventional acetabular cup systems.
- Increased joint stability, compared to conventional acetabular cup systems.

10. Implant Materials

Please refer to the relevant surgical technique associated with the system and the identification on the packaging for further information on implant materials. Further information on the material compositions is available from the manufacturer upon request.

- Cobalt-based alloy, CoCrMo, ISO 5832-12; ASTM F1537, (EndoDur-S)

11. Instrument Materials

Please refer to the REF list in this document and to the relevant surgical technique associated with the system and the identification on the packaging for further information on instrument materials.

Further information on the material compositions is available from the manufacturer upon request.

12. CMR Substances

Some system components contain cobalt as an alloy ingredient in a concentration above 0.1 % weight by weight.

Cobalt is listed as a CMR (carcinogenic, mutagenic and toxic to reproduction) substance.

The hazard class and category code(s) for cobalt are:

- Carc. 1B
- Repr. 1B

For identification of the affected components and further material specification, please refer to the REF list and the legend of materials in this document.

13. Implant Selection, permissible Combinations

For definitive identification information on the product such as system compatibility, article number, material and shelf life, refer to the identification on the implant and / or the packaging.

Also refer to the following sections of this document and to the relevant surgical technique associated with the system for further information on implant selection and permissible combination and for information on allocating and handling the instruments to be used for the implantation.

Combinations with implants from other manufacturers and / or combinations with LINK implants that deviate from the surgical technique specifications have not been tested and are not permitted.

14. Implant Diameters

Please refer to the relevant surgical technique associated with the system for further information on implant diameters.

15. Implant Anchoring

The Modular Trauma Head is implanted without cement.

16. Implant Lifetime

The lifetime of our implants is limited in principle and is determined by individual factors such as, for example, body weight and the level of activity of the patient, as well as by the quality and professional execution of the implantation. Based on these individual influencing factors, Waldemar Link defines the overall average lifetime of an implant based on its survival rate (i.e. the proportion of functional implants after a certain period of time starting from the time of implantation). According to the results of the tests performed, the survival rate of our implants corresponds to the general state of the art at the time of approval of the implants.

17. Reprocessing / Reuse

The implants are supplied as sterile single-use devices. Implants whose protective packaging is opened or damaged, or implants which have already been implanted are not permitted to be reprocessed or reused.

The following risks may occur if implants are re-used:

- Infections
- Reduced implant lifetime
- Increased wear and wear debris complications
- Disease transmission
- Inadequate implant fixation
- Limited implant function
- Implant response and / or rejection

Instruments must be disinfected and sterilized prior to use. For more information please refer to the description in the reprocessing instructions H50.

Additionally observe our separate packaging and cleaning instructions for instruments.

Single-use products may not be reused.

18. Resterilization

The Implants are designed for single-use only. Resterilization is not permitted.

Implants, as well as their materials are not suitable to be resterilized.

Unpredictable degradations may occur in these implants during resterilization.

For sterilization of the instruments, please refer to the description in the reprocessing instructions H50.

19. Storage and Transportation of Implants

Sterile-packaged implants must be stored in the undamaged original packaging in buildings with adequate protection against damage due to impacts, frost, humidity, excessive heat, and direct sunlight.

Further information is available from the manufacturer upon request.

20. Information for Patient Advisory

If the implantation of this system is considered to be the best solution for the patient and one of the circumstances described in section 21 is applicable to the patient, it is necessary to advise the patient with regard to the anticipated effects that these circumstances could have on the success of the operation. It is further recommended that the patient be informed about measures that he or she can take to reduce the effects of such complications. All information provided to the patient should be documented in writing by the operating surgeon. An implant ID for the patient must be handed over by the surgeon or hospital and the patient must be informed about the availability of a special patient information.

The patients should also be instructed:

- in detail about the surgery-related risks.
- in detail about the limitations of the implants, especially about the effects of excessive stress caused by body weight and physical activity, among other things. They should be encouraged to adjust their activities accordingly.
- about possible postoperative complications.
- about the material composition of the implant.
- that implants may respond to metal detectors during security checks (e.g. at airports) and carrying an implant ID as a proof is recommended.
- that implants may interact with medical imaging technique (e.g. MRI)

21. Circumstances that can interfere with the Success of an Operation

- Severe osteoporosis
- Severe deformities
- Local bone tumors
- Systemic diseases

- Metabolic disorders
- Case history of infections and falls
- Drug dependency or abuse, including excessive alcohol and nicotine consumption
- Obesity
- Mental disorders or neuromuscular diseases
- Heavy physical activities associated with strong vibrations
- Hypersensitivities

22. Warnings / Precautions

- The reuse of LINK single-use products is not permissible.
- Implants must be handled with great care and should not be modified or changed, even the smallest scratches and damages can considerably impair their stability or performance. Damaged implants are not permitted to be used.
- Surfaces provided for the connection of modular prosthetic components (cone, pins, screws) must not be damaged and may need to be cleaned with sterile liquid and dried before being joined together, so that neither blood nor any other coating impairs any of the connections, which could compromise the reliability of the connection.
- Do not manipulate or misuse instruments. We do not accept liability for products that have been modified, subjected to unintended use, or used improperly.
- For the processing of LINK instruments, it is presumed that the personnel have technical knowledge level I (Germany) and in other countries technical knowledge and expertise.
- Medical devices that are sent in for servicing must be processed beforehand in such a way that they cannot constitute a hazard to third parties.

23. Preoperative Planning

Preoperative planning provides important information to identify the appropriate implant system and select the components of a system. Make sure that all components required for the operation are laid out and ready in the operating room. Trial implants to verify proper fit (where applicable) and additional implants should be kept at the ready, in case other sizes are needed or the intended implant cannot be used. All LINK instruments necessary for the implantation must be on-hand, sterilized and intact.

If prosthesis implantation is indicated, then it must be taken into consideration, along with the overall circumstances of the patient:

- that all non-surgical and surgical treatment alternatives for the joint disease have been considered
- that artificial joint replacement performance is categorically inferior to natural joint performance, and an indication-related improvement in the preoperative condition is the only aim here
- that proper selection, placement and fixation of the devices are decisive factors, which will determine the life of the implant.
- that an artificial joint may loosen due to stress, wear and tear, and infection, or luxation or dislocation may occur
- that revision surgery, which under certain circumstances may exclude the possibility of restoring joint function, may be necessary due to loosening of the implant
- that the patient consents to undergo the operation and accepts the risks involved
- that if load-transferring bone structures are damaged, then the loosening of the components, bone and implant fractures, as well as other serious complications cannot be ruled out
- that if the patient is suspected of having allergies and tests positive on the applicable tests, then the patient's foreign body sensitivities (material tolerances) must be examined
- that acute and chronic infections, local and systemic, may compromise the successful implantation, so pre-operative microbiological analysis is recommended

Generally, the mechanical failure or fracture of an implant is a rare exception. However, this cannot be excluded with absolute certainty despite the sound structure of the implant.

This may be due to stress on the implant and prosthesis as the result of a fall or accident, among other things.

If the bone area where the implant is anchored is altered in such a way that the prosthesis is no longer able to withstand normal stress and an area of the prosthesis becomes subject to a stress imbalance, then a mechanical failure of the implant system may result. Such stress imbalances may also occur if the anchoring elements of the implants are required to form a bridge over larger bone deficiencies without optimal reinforcement of the bone. It is recommended that the implant with the largest possible anchoring elements be used.

Proper preparation for surgical procedures also includes the functional testing of implants and instruments prior to use.

For definitive identification information on the product such as system compatibility, article number, material and shelf life, refer to the identification on the implant and / or the packaging. You should also take advantage of the training courses and printed materials provided for your information. To learn more, please contact the Waldemar Link GmbH & Co. KG sales office or your field representative.

24. Handling

All implant components are supplied sterile as single-use devices in individual packages. The implant components are sterilized by gamma sterilization, at least 25 kGy. By contrast, components made of highly cross-linked polyethylene or highly cross-linked polyethylene with vitamin E are sterilized with ethylene oxide (ETO).

Implants should always be stored in their unopened protective packaging. Examine the packaging for damage before using the implant. Damaged packaging can have an adverse effect on both the sterility of the device as well as the proper performance of the implant, such that the device may no longer be used.

- Check the use by date on the implants. Implants with expired use by dates are no longer permitted to be used for implantation!
- After opening the package, check that model and size of the implant are matching with the information printed on the package labeling.
- Observe the pertinent standards for the aseptic handling of devices during and after removal of the implant from the packaging.
- When removing the packaging, make a record of the batch or serial numbers on the label, since this information is decisive for batch tracing. Self-adhesive labels with this information are enclosed with every package for your convenience.

Instruments must always be treated with care, this particularly applies during transport, cleaning, maintenance, sterilization, and storage. The sterile status of the instruments depends, inter alia, on the sterile items packaging and the prevailing storage conditions and must be established together with the operator's hygiene officer on a case-by-case basis. Direct sunlight must be avoided. Improper handling and care, as well as unintended use, can lead to premature wear or damage.

Devices made of plastics (e.g. PP-H) may not be located by means of an external imaging device.

25. Intraoperative Use

Please refer to the relevant surgical technique associated with the system for information about the intraoperative use of the system.

26. Postoperative

In addition to movement and muscle training, special attention must be paid to carefully instructing the patient during the postoperative phase.

Physician-supervised postoperative monitoring of healing progress is recommended. Where applicable, patients should also be advised on how to avoid overstraining themselves.

Follow-up examinations should be carried out regularly or immediately if symptoms occur.

27. Notes on MRI and CT Examination Procedures

Our implants were not evaluated for safety and compatibility for MRI and CT examination procedures.

In the case of our metallic implants and implant components, MRI examinations pose potential risks to the patient due to possible heating and migration of the implants or implant components.

Likewise, there is a potential risk of artifact formation in MRI and CT examinations of our metallic implants and implant components.

The probability of occurrence and the extent of the potential risks mentioned depend on the type of device used, its device parameters and the sequences used.

Always follow the instructions in the operating instructions of the manufacturer of the device used for the imaging.

The selection of the imaging examination procedure and the assessment of possible side effects is the responsibility of the examining physician.

The examining physician must take into account the individual condition of the patient and other diagnostic methods.

28. Explantation of Implants / Revision Surgery

Please refer to the relevant surgical technique associated with the system for information about explantation of the implants and revision surgery.

29. Disposal

Packaging and system components to be discarded must be handled in compliance with your national and local regulations for hospital disposal.

30. Instruments

Please refer to the description in the reprocessing instructions H50 for the:

- Initial use
- Performance test
- Maintenance
- Manual cleaning
- Cleaning in a washer-disinfector
- Sterilization
- Servicing

31. Link to Summary of Safety and Clinical Performance

At the time of creation of this document, the EUDAMED database was not yet active, therefore no link to the summary of safety and clinical performance can be given here.

32. Requests

Requests of any kind should be directed to Waldemar Link GmbH & Co. KG (see contact information in this document).

33. Complaints about our Products

All complaints must be addressed to Waldemar Link GmbH & Co. KG at: complaint@link-ortho.com.

In the event of a complaint, the name or reference number of the corresponding component should be specified with the serial number (SN) or the lot number (LOT), your name, and your contact address. The reason for the complaint should be given in brief.

34. Report of serious incidents

Any serious incident that occurs in relation to the device must be reported to the manufacturer and the authority responsible for your location.

Waldemar Link GmbH & Co. KG and/or other corporate affiliated entities own, use or have applied for the following trademarks in many jurisdictions: LINK, EndoDur

Other trademarks and trade names may be used in this document to refer to either the entities claiming the marks and/or names or their products and are the property of their respective owners.

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